


Infant Test Lung

Neonate and infant test lung




Test your neonatal or infant ventilator with the **imtmedical** Infant Test Lung.

Features

- Two different compliances
- Two different resistances
- Two separate test lungs
- Individual lung selector
- Leak simulator on each lung
- Simple to use

Specification

- **Tidal Volume:**
150ml or 300ml
- **Dimensions:**
120mm x 250mm x 60mm
- **Compliance:** constant over the full range at 0.35 ml/mbar or 0.70 ml/mbar:



Ordering Information

Partnumber: 500.040.000

For further information contact:

imtmedical ag
Gewerbestrasse 8
9470 Buchs – Switzerland

Tel: +41 81 750 66 99

Fax: +41 81 750 66 95

Mail: sales@imtmedical.com

Internet: www.imtmedical.com